

# How To Make TE Buffer pH 8.0

## About TE buffer

TE buffer (Tris-EDTA) is a commonly used buffer solution for resuspending and storing nucleic acids, especially DNA. The Tris solution keeps the DNA soluble in water while EDTA, a chelator of cations such as magnesium, protects nucleic acids against enzymatic degradation.

## Recipe

The recipe below is to create a 100 mL TE buffer solution.

Reagent	Volume	Final concentration
1M Tris-Cl (pH 8.0)	1 mL	10 mM
0.5M EDTA (pH 8.0)	0.2 mL	1 mM
Distilled H <sub>2</sub> O	98.8 mL	

## How to make TE buffer

1. Measure out 1 mL [1M Tris-Cl \(pH 8.0\)](#) and add to a 100 mL Duran bottle.
2. Measure out 0.2 mL [0.5M EDTA \(pH 8.0\)](#) and add to the Duran bottle.
3. Top up the solution to 100 mL by adding 98.8 mL of distilled water.
4. Place the lid on the bottle and invert a few times to mix.
5. To sterilise, autoclave the solution on a liquid cycle (20 min at 15 psi).

## Storage of TE buffer

Store TE buffer at room temperature (+15°C – +25°C).

## Safety

TE buffer is not classified as hazardous. However, always be sure to read the TE buffer [safety data sheet](#) before use.

### Steven Bradburn, PhD

Steven is the founder of Top Tip Bio. He is currently a Medical Writer and a former Postdoctoral Research Associate. Enjoyed the tutorial? Then let me know by leaving a comment below, or consider [buying me a coffee](#).

